

Advantages of the HEARING AID SWEAT BAND™ over other BTE hearing aid moisture protectors:

- Not sprayed inside the hearing aid. There is nothing that will harm the inner controls.
- You'll be able to operate the controls with the HEARING AID SWEAT BAND™ in place.
- Has the ability to absorb large amounts of moisture enabling the wearer to leave in their hearing aid during a visit to the hair salon. In addition, it protects the BTE from hair spray and other hairstyling products.
- Easy to slip over the hearing aid; be sure to order it long enough to cover the microphone. No need for an application procedure or special tools.

- Protects the BTE from wind and sand at the beach, wind and snow on the slopes, wind and water on the boat.

- Non-irritating and soft to the skin.



"Hearing aid repair technicians have found that one of the most common causes of hearing aid breakdown is moisture. What I found most interesting about the HEARING AID SWEAT BAND™ was its ability to hold a large amount of moisture."

~Michael J. Phillips,
Education Coordinator Lori/Unitron

Why use the HEARING AID SWEAT BAND™?

Because...

- Made from specially designed fabric, the HEARING AID SWEAT BAND™ protects the BTE from the damaging effects of moisture and dirt infiltration.
- The special fabric is hypoallergenic, soft and supple and is not irritable to the skin.
- The HEARING AID SWEAT BAND™ allows sound to readily pass through it. Tests have shown the sound dampening effect to be negligible.
- Allows moisture to evaporate and wick away. Doesn't hold moisture against the hearing aid.
- The HEARING AID SWEAT BAND™ reduces cleaning and maintenance visits.
- Reduces the damaging effects of moisture on the batteries.
- The HEARING AID SWEAT BAND™ helps protect the inner workings of the analog and digital BTE.



Hearing Aid
SWEAT BAND™



1-866-644-2500

Fax: 1-888-798-2263

www.hearingaidsweatband.com

Copyright © 2007 VanB Enterprises

DOES YOUR BTE FAIL DUE TO PERSPIRATION, MOISTURE, WIND OR DUST?



THE HEARING AID SWEAT BAND™ IS THE ANSWER.

VanB Enterprises
www.hearingaidsweatband.com

THE HEARING AID SWEAT BAND™

STANDARD WIDTH	
Mini	1"
Small	1¼"
Medium	1½"
Large	1¾"
X Large	2⅞"
5 XL	3"
SLIM SERIES	
Mini Slim	1"
Small Slim	1¼"
Medium Slim	1½"

"We started recommending the (HEARING AID SWEAT BANDS™) to our clients...they do a great job!"

~Betty L. Jensen
M.A., CCC-A,
a certified clinical audiologist
and owner
of Tri-State
Hearing Centers



Virtually anyone who wears a BTE can benefit from a **HEARING AID SWEAT BAND™**.

Great for runners, golfers, walkers, skiers, soccer and tennis players, helmet wearers, children, and anyone who works outdoors or in a dusty environment.


The **HEARING AID SWEAT BAND™** eliminates sudden loss of hearing due to moisture from perspiration, rain or snow, and stops dust infiltration. The **HEARING AID SWEAT BAND™** also offers advanced improvement over feedback due to wind.

Controls can be easily operated even with the **HEARING AID SWEAT BAND™** in place. The protective covering simply slips over the hearing aid, no need for an application or special tools.



Slide fully over BTE to cover microphone area

Hearing aid repair technicians have found that one of the most common causes of hearing aid breakdown is moisture and dust infiltration. The **HEARING AID SWEAT BAND™** offers the solution.



The **HEARING AID SWEAT BAND™** is made from a specially designed fabric that is both moisture repellent and moisture absorbent but does not inhibit sound. It's inexpensive and disposable, but will withstand a few launderings. We recommend hand washing with a mild detergent and air drying.

The **HEARING AID SWEAT BANDS™** are conveniently packaged with two per pack and are available in six colors and nine sizes ranging from Mini to 5 XL and now the new "Slim Series."

The "Slim Series" of the **HEARING AID SWEAT BAND™** keeps pace with the new smaller Behind the Ear hearing aids on today's market. The "Slim Series" is available in Mini, Small and Medium lengths in Natural, Light Brown and Gray; and provides the same protective benefits as the standard width **HEARING AID SWEAT BANDS™**.

The **HEARING AID SWEAT BAND™** is also available in 5 XL which was designed to protect FM Systems allowing FM System wearers to now have the capability of keeping their system safe from moisture, perspiration, dust and other external contaminants. Designed to meet the specific needs of FM Systems, the 5 XL measures 3" in length and comes in all colors.

